

Specifications

Shimmer3 ECG/EMG Unit



User Guides
Sample Data
Case Studies
Software
Expansion Boards

Key Features

- 5-wire, 4-channel ECG, measuring bipolar limb leads and user's choice of V1 - V6
- Measure 2 channels of EMG data with a common reference electrode
- Software configurable right-leg drive for common-mode interference rejection
- Software configurable amplifier gain
- Software configurable data rate
- Respiration demodulation capability on-chip
- Lead-off detection capability on-chip
- Test signal on-chip for validation purposes
- EEPROM storage device (on the ExG daughterboard) enables expansion board detection and identification, as well as 2032 bytes of data storage available to user

Introduction

The ECG/EMG unit, previously also known as the ExG unit, provides a configurable digital front-end, optimised for the measurement of physiological signals, for example 5-wire ECG (Electrocardiography) and 2-channel EMG (Electromyography).

Compatible with the Shimmer3 platform, the ECG/EMG unit also boasts best data quality with integrated 10 DoF inertial sensing via accelerometer, gyroscope, magnetometer and altimeter, each with selectable range.

The ECG/EMG unit also provides highly accurate and scientifically reliable raw data to allow complete control over capture and interpretation of sensed data in real-time.

Product Overview

While addressing the challenges of mobility, the ECG unit can record the pathway of electrical impulses through the heart muscle, and can be recorded on resting and ambulatory subjects.

When configured for EMG, the unit can also measure and record the electrical activity associated with muscle contractions, assess nerve conduction, muscle response in injured tissue, activation level, or can be used to analyse and measure the biomechanics of human movement.

Applications

Atrial fibrillation
Premature ventricular contraction
Heart function monitoring
Abnormal rhythm detection and alert
Biomechanics, muscle activity, gait and posture disturbance
Fatigue analysis
Sports technique, performance and medicine
Neuro Rehabilitation
Tremor Analysis
Veterinary Science
Orthopaedics
Thoracic Bioimpedance (Respiration)

Technical Specifications

Gain	Software configurable (1, 2, 3, 4, 6, 8, 12)
Data Rate	Software configurable (125, 250, 500, 1000, 2000, 4000, 8000 SPS)
Input Differential Dynamic Range	Approx 800 mV (for gain = 6)
Bandwidth	8.4 kHz
Ground	Wilson Type Driven Ground
Input Protections	ESD and RF/EMI filtering; Current limiting; inputs include defibrillation protection (survive only, not repeat)
Connections	EMG: Input Ch1N, Input Ch1P, Input Ch2N, Input Ch2P, Reference (Ref) ECG: Input RA, Input LA, Input LL, Input Vx, Reference (RL) All Hospital-Grade 1mm Touchproof IEC/EN 60601-1 DIN42-802 jacks
Weight	31 grams
Dimensions	65 x 32 x 12 mm
EEPROM Memory	2048 bytes

Shimmer3 Specifications

Processing	MSP 430 microcontroller (24mHz, 16Bit)
Communication	Bluetooth – RN4678
Storage	Integrated 8GB microSD card slot
Battery	450mAh rechargeable Li-ion
Integrated 3-Axis Accel.	ICM-20948
Accel. Range	±2g, ±4g, ±8g, ±16g

ECG Module 5-wire, 4- channel ECG, measuring bipolar limb leads and user's choice of V1 - V6
Digital interface includes test signal generation for validation purposes
Respiration demodulation from ECG data and lead-off detection

EMG Module Two channels of EMG data
Digital interface includes test signal generation for validation purposes

Supporting Software Shimmer ConsensysPRO & ConsensysBASIC
Shimmer LabVIEW Instrument Driver
Shimmer MATLAB Instrument Driver
Shimmer Java / Android API
Shimmer C# API
Shimmer 9DoF Calibration

Supporting Hardware & Accessories Biophysical Leads
Straps, Documents, Charging Dock/Base, Case
Electrodes
Shimmer3 Calibration Stand
